PULSOSSIMETRO PORTATILE

SANICO

PL101

IT EN FR ES

PULSOSSIMETRO PORTATILE

IT

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FINGERTIP PULSE OXIMETER

EN

This digital device is a non-invasive, compact, practical medical device with low energy consumption. Just insert the finger tip into the sensor of the device and you will read the SpO2 values, the pulse frequency, the pletismographic wave (display 1-2) and the perfusion index (ratio between pulsatile blood flow and static flow in peripheral tissues such as the fingertip; it indicates the strength of the pulse to the sensor: values 0÷20).


GENERAL DESCRIPTION

Haemoglobin Saturation is the percentage between the capacity of Oxyhemoglobin (HbO2) that compounded with oxygen and that of all combined hemoglobin (Hb) obin (HbO2) in blood. In other words, it is the saturation of Oxyhemoglobin in blood. It is a very important physiological parameter for Respiratory and Circulation Systems. Many respiratory diseases could reduce haemoglobin saturation in human blood. Moreover, factors such as Automatic Organic Regulation Malfunction caused by anaesthesia, trauma resulted from major operation and some medical examination can also cause problems in oxygen supply, which might reduce human haemoglobin saturation. As a result, such symptoms as megrim, vomiting and asthenia might appear to patients. Hence, it is very important to know hemoglobin saturation of patient timely in clinical medical aspects. The fingertip pulse oximeter features in small volume, low power consumption, convenient operation and portability. It is only necessary for patient to put one finger into fingertip photoelectric sensor for diagnosis, and the display screen will directly show measured value of hemoglobin saturation. It has been proved in clinical experiments that it possesses rather high precision and repeatability.

MEASUREMENT PRINCIPLE

The principle of the oximeter is as follows: An experience formula of data process is established by exerting Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin(RbH) and oxyhemoglobin (O2 Hb) in glow and near-infrared zones. Operation principle of the instrument is to combine Photoelectric Oxyhemoglobin Inspection Technology with Capacity Pulse Scanning and Recording Technology, so that two lights with different wavelength (660nm glow and 940nm near infrared light) can be focused onto human nail through perspective clamp finger-type sensor, and then measurement signal got via a photosensitive element. Relevant measurement information will be processed in electronic circuits and microprocessor and then showed up on the OLED display.

- Generating probe

- Receiver photocell

PRECAUTIONS FOR USE

- The device does not have an alarm indicating abnormal values. Do not use it in situations where alarms are required. Do not use the pulse oximeter in an MRI or axial tomography environment.
- Explosion hazard: do not use the pulse oximeter in explosive atmospheres.
- The use of the pulse oximeter is indicated only as an aid in patient evaluation. The device should be used in conjunction with other methods for the evaluation of clinical signs and symptoms.
- Frequently check the area of application of the pulse oximeter to determine the positioning of the sensor and the correct circulation and sensitivity of the patient's skin.
- The finger pulse oximeter is indicated for both home use and hospital use.
- The SpO2 sensor and the photoelectric reader must be positioned so that the patient's artery is in between them.
- Prolonged use or patient status may require you to change the sensor

application site regularly. Vary the sensor application site and check skin integrity, circulatory status and correct orientation at least every 4 hours.

- Inaccurate measurements could be caused by autoclaving, ethylene oxide sterilization, immersion of the device in liquids, significant levels of dysfunctional hemoglobins (such as carboxi-hemoglobin or methemoglobin), intravascular dyes such as indocyanine green or methylene blue.
- SpO2 measurements could be negatively affected in the presence of high ambient light (avoid direct exposure of the sensor to sunlight), excessive patient movements, vein pulsations, placement of the oximeter on a finger and simultaneous presence of a pressure cuff on the arm, from the presence of an arterial catheter or an intravascular line, in patients with hypotension, severe vasoconstriction, with severe anemia or hypothermia, in the presence of cardiac arrest or in a state of shock.
- The presence of nail polish or false nails can cause inaccurate readings.

FEATURES

- OLED display
- Product adopts double color OLED display can show the six different display mode
- Low-power consumption, continuously work for 40 hours
- Low voltage indicator
- In the absence of signals, the product will be in after 10 seconds to enter standby state
- Visual and auditory prompt function.
- Small in volume, light in weight, and convenient to carry
- Product Operation Scope
- The fingertip oximeter can be used to measure human haemoglobin saturation and heart rate through finger. The product applies to using in family, hospital (including clinical use in internist/surgery, anaesthesia, paediatrics, intensive care, etc.), oxygen club, social medical organizations, physical care in sports (It can be used before or after sports. Operation in sport procedure is not recommended). It is also applicable to enthusiasts on mountaineering, patients (convalescents at home or those need first aid treatment), elders over 60, those work more than 12 hours, sports and those work in the hermetic circumstance, etc. The product is not suitable to monitor patient continuously.

OPERATION INSTRUCTIONS


- Hold the lock button on the battery cover, push the battery cover down at the same time, the battery cover is opened.
- Install two AAA batteries into battery cassette before covering its cover.
- Plug one finger into rubber hole of the oximeter (it is best to plug the finger thoroughly) before releasing the clamp with the nail upwards.
- Press button on the front panel;
- Don't tremble your finger when the oximeter is working. It is recommended to stay standstill while using your oximeter.
- Press the button on the front panel, if we want change display direction;
- Read relevant datum from display screen.
- The instrument has the function of sleep, no signal will enter standby state of sleep;
- Please replace new batteries when OLED indicates the batteries are in low power.

When plugging your finger into the oximeter, your nail surface must be upward.

Note: Please use medical alcohol to clean the rubber before each test and clean the tested finger with alcohol before and after the test. (The rubber inside of the oximeter adopts medical rubber, which has no toxin, no harm, and brings no side effect such as allergy to the skin).

BRIEF DESCRIPTION OF FRONT PANEL

ON / OFF button and display options

When the measurement is displayed, pressing the ON / OFF key changes the orientation of the display to make reading easier: in positions 1 and 2 the Pletismographic Wave (pulse wave) is also visible.

DETAILED DESCRIPTIONS OF PRODUCT FUNCTIONS:

- Display type: OLED
- SpO2 measurement range: 70% -99% Precision: for SpO2 values between 70% -99% $\pm 2\%$; for values ($\leq 70\%$) not specified Resolution: $\pm 1\%$
- Measurement range PR: 30BPM-240BPM Accuracy: $\pm 1\text{BPM}$ or $\pm 1\%$ (greater)
- Power supply: two 1.5V alkaline batteries (AAA)
- Supply current: <30mA
- Automatic shutdown: after 10 seconds without signal.
- Environmental characteristics: Working temperature: 5° C-40°C Storage temperature: -10° C-40°C Working humidity: 15% -80% Environmental conservation humidity: 10% -80% Working pressure: 70kPa -106kPa

The device complies with IEC 60601-1-2 for EM emissions. The materials which the user can come into contact with have no toxicity and no action on the tissues and are in accordance with ISO 10993-1, ISO 10993-5 and ISO 10993-10.

Subject to technical modifications.

CLASSIFICATION

Management Class for Medical Devices: equipment Anti-electric Shock Type: Internally powered equipment Anti-electric Shock Degree: Type BF equipment

MAINTENANCE AND PRESERVATION

- Replace the batteries timely when low voltage lamp is on.
- Clean the surface of fingertip oximeter before it is used to diagnose patients.
- Remove the batteries inside if you will not operate the oximeter for a long time.
- It would be better to preserve the product in -10-40 (14-104) and humidity is 10%-80%.
- It is recommended that the product should be kept dry anytime. A wet ambience might affect its lifetime and even damage the product.
- Please follow the law of the local government to deal with used batteries; subjected to WEEE Directive, dispose of according to local laws.

PRODUCT ACCESSORIES

- Lanyard and bag for carrying
- Two 1.5V AAA alkaline batteries
- User Manual

GUIDANCE AND MANUFACTURE'S DECLARATION - ELECTROMAGNETIC RADIATION - FOR OTHER EQUIPMENTS AND SYSTEMS

The Pulse Oximeter is designed to be used in specified electromagnetic environment . Users of the Pulse Oximeter must use it in the following environments.

Radiation Test	Compliance	Electromagnetic environment-guidance
RF interference	Group 1	RF signal of Pulse Oximeter is simply created by its internal function. Therefore, its RF interference is very low and is not likely to cause any interference to nearby electronic equipment.
RF interference CISPR 11	Class B	The Pulse Oximeter applies to all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

POSSIBLE PROBLEMS AND RESOLUTIONS

Problem	Possible reason	Solution
SpO2 or PR can not be shown normally	1. Finger is not plugged correctly 2.Patient's Oxyhemoglobin value is too low to be measured	1. Retry by plugging the finger 2. Try more times. If you can make sure there is no problem in the product, please go to hospital timely for exact diagnosis
SpO2 or PR is shown unsteady	1. The finger might not be plugged deep enough 2. Finger is trembling or the patient is on movement status	1. Retry by plugging the finger 2. Please remain at rest
The oximeter cannot be turned on	1. Inadequate power or power off 2. Batteries might be installed incorrectly 3. The oximeter might be damaged	1. Replace the batteries 2.Please reinstall the batteries 3.Please contact with local customer service centre
Indication lamps are suddenly off	1. The product automatically shuts off when no signal is detected in 10 seconds 2. Inadequate power	1. Normal 2. Replace the batteries

SYMBOLS AND DEFINITIONS


Symbol	Definition
	The equipment type is BF
	Refer to user manual before use
%SpO2	Hemoglobin Saturation
PR BPM	Heart Rate(BPM)
	Low power indication
	Batteries and electronic devices must be disposed of in accordance with local regulations
SN LOT	Serial Number/ Lot Number
EC REP	Prolinx GmbH - Germany
	Xuzhou Yongkang Electronics Science Technology Co., Ltd - China
CE 0123	Medical Device certified in accordance with Dir. 93/42/CE Class IIb

HOW TO REPLACE BATTERIES:

- To remove the back cover compartment push the button and follow the printed arrows direction
- Install two AAA batteries into the battery compartment. Match the plus (+) and minus (-) signs in the compartment. If the polarities are not matched, damage may be caused to the oximeter.
- Slide the battery door cover horizontally along the arrow shown as the picture.

Notes:

- Please remove the batteries if the pulse oximeter will not be used for long periods of time.
- Please replace the batteries when the power indicator starts flickering.

USING THE LANYARD

- Insert the thinner end of the cord through the hole in the device.
 - Pass the other end of the cord through the eyelet, then pull.
- Warnings!
- Keep the device away from children. Some parts are small and could be swallowed.
 - Pay attention to the risk of strangulation in the presence of cables or pipes.
 - The cord tied to the device can cause strangulation due to excessive length.

WARRANTY

This appliance is covered by a 3-year (2 years European warranty plus 1 additional year) warranty from the date of purchase for defects in material or workmanship. Warranty is valid only upon presentation of sales receipt or invoice. Opening or tampering will void the warranty. Accidental damage caused by improper use or non-compliance with the instructions is not covered. For info contact your local distributor.

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Xuzhou Yongkang Electronic Science Technology Co., Ltd
Company name: Prolinx GmbH
Imported by: SANICO Srl - Via G.Ferraris, 37/41