SANICO

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INSTRUCTION MANUAL

Please read this instruction manual carefully
before using the device

The manufacturer declares that the product complies with the European Union's Medical Device Regulation (EU MDR) 2017/745 Class IIa.

Contents

Introduction

Intended for Use

Clinical benefit

User restrictions/Contraindications

Important Information before Use

Product Identification

Instruction for Use

Instruction for Use (Child users)

Battery Recharge

Change Air Filter

Cleaning and Disinfecting

Technical Specification

Consult your physician or doctor if

If you experience these side effects

Seek immediate medical attention

Troubleshooting

EMC Tables

Explanation of Symbols

Introduction

This product is intended to be used for inhaling medication like aerosol treatment for respiratory disorders such as mild asthma.

- 1.Easy to operate.
- 2.Compact and travel convenient
- 3.Allow to use with face mask or mouth piece
- 4.Nebulization rate ensure optimized treatment duration
- 5.Rechargeable battery operated: more eco-friendly
- 6.Particle size of 2–5 μm are optimal for reaching around the lungs
- 7.This device is intended for use by patients in homecare and healthcare settings.
- 8.Ask your doctor or physician for further information about the potential applications.
- 9.The device does not need to be calibrated.
- 10.The device does not supply with medication.

Intended User

This product is a portable aerosol nebulizer device that operates electronically with an internal rechargeable lithium polymer battery.

The device functions by creating compressed air using its internal pump, and this compressed air flows into the nebulizer kit. Within the nebulizer kit, the air pressure forces the liquid medication through the vaporizer head, where it is converted into an aerosol.

This aerosolized medication is then delivered for inhalation through either a mouthpiece or a mask, providing an effective means of treatment.

Operating Principle

The operating principle of the device involves pumping the medication through a dedicated channel, where it is mixed with the compressed air generated by a compressor pump. This combination of medication and compressed air results in the transformation of the medication into fine particles, which are then sprayed into the air when they come into contact with the compressed air. This process ensures that the medication is in an inhalable form for effective delivery to the patient's respiratory system.

Intended for Use

This product is designed delivering prescribed medications to manage respiratory conditions like asthma in both the homecare and healthcare settings. The nebulizer kit is intended for multiple uses by single patient.

The device is compatible and suitable for delivering these specific common asthma medications such as Albuterol/ Salbutamol and Terbutaline. For other prescriptions meds please obtain advice from your doctors or physician to ensure safe and effective use of the nebulizer with any other prescribed drugs.

Indication for Use

The nebulizer provides significant relief to individuals by delivering medication directly to the lungs, reducing airway inflammation, and improving lung function.

Clinical benefit

The intended clinical benefits of the Sanico nebulizers include:

- Reduce drug dosages
- Accelerate onset of action
- Reduce side effects (such as airway inflammation)
- Enhance convenience in comparison to systemic drug delivery

Intended User

1. Users use nebulizers on the advice of doctor to provide relief from asthma symptoms, range from age 2 and above up to 65 years old.
2. For children who are unable to perform the aerosol treatment by themselves, require assistance from their guardian or parents. The guardian or parents should fully understand the instruction manual, including how to operate the nebulizer and all precautions and warnings associated with its use.
3. Both professional healthcare personnel and laypersons.

User restrictions/Contraindications

If a patient has any of these contraindications, alternative treatment options or adjustments to their treatment plan may be necessary.

- Unconscious patients or patient not breathing spontaneously.
- Heart rhythm disorder, cardiac irritability (may result to dysrhythmias)
- Internal bleeding
- Patients with unstable and increased blood pressure/pulse
- Individual intolerance or allergic reaction to Albuterol/ Salbutamol and Terbutaline other medications, foods, dyes, or preservatives
- Diabetes (high blood sugar)
- Seizures
- Thyroid disease
- Pregnant or trying to get pregnant, miscarriage
- Breast-feeding
- Liver disease
- Osteoporosis, weak bones
- Immune system problems
- Uterine infection
- Eye disease, such as glaucoma, cataracts, or blurred vision

WARNING:

- If you take too much dose, seek emergency medical attention straight away even if you do not feel unwell or have symptoms of an overdose. An overdose can be fatal. Contact your doctor immediately or go to the Emergency Department of your nearest hospital.
- If you take this medication on a regular basis, take it as soon as you can. If it is almost time for your next dose, take only that dose. Do not take double or extra doses.
- Aerosol treatment with drugs may affect your ability to drive or use machines. If you feel sick, dizzy or tired, or give you a headache or blurred vision, do not drive or use machines and contact your doctor immediately.
- If the user uses wrong medication or unknown medication, contact your doctor immediately or go to the Emergency Department of your nearest hospital .
- Do not mix your medication with other drugs in the same nebulizer. Consult your doctor if you also on other type of medication. Clean and disinfect the nebulizer and parts thoroughly if use other type of medication.
- The device does not act as a MDI medication dose counter and does not include a dose counting function. *MDI= Metered Dose Inhaler
- The device is not intended to diagnose, measure, cure medical condition or replace an interpretation from a health care provider; Patients are expected to always follow their doctor's advice.

Important Information before Use

1. Do not expose the device to extreme temperature, very high humidity level.
2. Avoid heavy object place on top of device, extreme stress, shock or accidental drop down.
3. You shall apply force on the power key, push upward to start the device or downward to turn off the device.
4. The nebulizer is for inhaling for respiratory disorders only. Do not use this device other than the intended purposes.
5. During the use of the device, the user should not hold his/ her breath.
6. If you feel anything unusual when using the device, stop and consult a doctor immediately.
7. Always dispose of remaining medication after use, and use fresh medication each time.
8. Store the parts in a clean location to avoid infection. Do not use or store the device where it may be exposed to noxious fumes or volatile substances.
9. Do not disassemble or attempt to repair the nebulizer or charging port.
10. If you would require any assistance during setting up, using or maintaining, please contact the local distributor.
11. The device is not suitable for use in an anaesthetic breathing system or a ventilator breathing system.
12. Do not use tap or mineral water in the nebulizer for nebulizing purposes.
13. During aerosol process, avoid turning off the power unexpectedly. Sudden power cut may results in treatment under dose, treatment interruption or treatment incomplete.
14. Do not use a damaged nebulizer kit, mouthpiece or facemask.
15. Make sure that the vaporizer head is correctly installed before use.

Important Information before Use

16. The nebulizer and charging cable line are not waterproof. Do not spill water or other liquids on these parts. If liquid does spill on these parts, immediately unplug the charging cable line and wipe off the liquid with gauze or other absorbent material.
17. To avoid the medication residue be sure to wipe the face after removing the mask.
18. Please keep this device away from sunlight.
19. The nebulizer may become hot during long period of operation.
20. Use only original nebulizing parts, air filter and air filter cover.
21. Do not add more than 5 ml of medication to the medication tank.
22. Do not carry or leave the nebulizer kit while the medication tank contains medication.
23. Do not use the device while sleeping or if drowsy.
24. Approved for human use only.
25. If any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user is established.
26. Recommend at least 2ml of prescription medication to operate the nebulizer.

WARNING:

1. You should always follow the instructions of your doctor regarding the type of medication to use, the dosage, and the frequency and duration of inhalation.
2. Only use medication prescribed or recommended by your doctor or physician.
3. Do not modify this equipment without authorization of the manufacturer.
4. Do not use the device where it may be exposed to flammable gas.
5. If you experience any allergic reactions or other difficulties during use, stop using the device immediately and consult your doctor.

Important Information before Use

1. Keep the device out of reach of unsupervised infants and children. The device may contain small pieces that can be swallowed.
2. Avoid potential allergic reactions; please avoid the front of tip in direct contact to patient's wound.
3. Always use the device in appropriate operating environment. If the device experienced extreme temperature/humidity/atmosphere, allow it to settle for 2 hours before use.
4. This product should not be used on patients, who are unconscious or are not breathing spontaneously. If you feel anything unusual when using the device, stop and consult a doctor immediately.
5. To completely shut down the device, remove the charging cable line from the charging port.
6. Never use the device near flammable gases, oxygen or nitrogen oxide
7. Do not leave the device unattended with infants, children or compromised individuals. Must have the supervision or guidance when using the nebulizer with infants, children for home treatment at all time.
8. Changes or modifications not approved by manufacturer will void the user warranty.
9. Make sure all the parts and accessory are not missing or broken, try to assemble and dismantle. If experiencing trouble assemble or dismantle, contact local distributor for detail instructions.
10. Do not cover the nebulizer with a blanket, towel, or any other type of cover during use. The nebulizer could get hot and there is a possible risk of a burn if touched.
11. When using this device, there will be some noise and vibration caused by the pump in the nebulizer. There will also be some noise caused by the emission of compressed air from the nebulizer kit. This is normal and does not indicate a malfunction.
12. Only fill the amount of medicine prescribed into the medication tank. Pay attention, that the maximum level will not be exceeded.
13. One prescribed medication at a time. Do not mix two or more kinds of prescribed medication at a time.
14. When use with mouthpiece, do not cover the exhalation port.

Product Identification

- 1. Mouthpiece
- 2. Inhalation Top
- 3. Medication Tank
- 4. Child Mask
- 5. Front Case
- 6. Adult Mask

- 7. Air Filter Cover
- *Air filter inside
- 8. Power Key
- 9. Vaporizer Head
- 10. Bottom Case

Instruction for Use

1. Inspect the device and parts before using them each time, and check the components make sure that there are no problems. In particular, be sure to check the following:

- That the Vaporizer Head is not damaged or blocked.
- That the nebulizer operates normally.
- The accessories are good to use.

2. Make sure that the parts are attachable and attached correctly.

3. Make sure that the air filter is correctly attached.


1. Take out the inhalation top part. Unscrew counter clockwise the inhalation top vaporizer head from the medication tank.

2. Add the correct amount of prescribed medication to the medication tank but do not add more than the MAX scale.


3. Once finished filling with prescribed medication, rotate clockwise to tighten the inhalation part to the medication tank. Make sure that the vaporizer head is securely attached to the inhalation top.

Instruction for Use

4. Attach the mask or mouthpiece to the inhalation top.

The mask is supplied by qualified supplier with verification of biocompatibility. Mask ID take in consideration with the nebulizer main body and inhalation top and mouth piece

5. Attach the inhalation part to nebulizer.

Clockwise to lock the device, counter clockwise to detach the inhalation top and the device.

6. Switch on the power button and make sure the nebulization begins to aerosol. While the nebulizer is operating, the user will have continuously inhalation and exhalation.

7. When the treatment is complete, switch off the power button.

For Child users, requires adults/guardian/caregivers supervision at all time.

For Child users who cannot understand the instruction manual, requires adults/guardian/caregivers to assist to perform aerosol treatment.


1. Please follow the previously mentioned assembly steps. After that, take out the child's mask, secure it onto the child's face, ensuring that the mask fits comfortably and covers both the mouth and nose. If the child feels uncomfortable, make necessary adjustments, and loosen the mask slightly.
2. Initiate the treatment and ensure that the device and child are under constant supervision. Do not leave the child unattended while using the device.
3. After completing the treatment, gently clean the child's face using a towel or wipe with clean water.
4. After the treatment, dispose of any unused medication as per disposal guidelines. Then, follow the cleaning and disinfection steps to clean the device immediately to maintain its hygiene.
5. Make sure to store the device and any medication in a safe and secure location where they cannot be accessed by children or pets.
6. Before using the device again, carefully examine it and ensure there are no medication residues. Also, make sure the device is completely dry before the next use to prevent any issues.

Battery Recharge

The device is built in with lithium rechargeable battery, do not dismantle or replace the battery with our manufacturer's authorization. The blue indicator lights up when charging. Turns green light when fully charged.

This device is powered with a built-in lithium battery. The battery can be charged with the USB cord. (Located back of the device)

1. Connect Type-C end of USB cord to the device.
2. Connect Type-A end of UCB cord to the adapter.
3. The charge indicator light shows "BLUE" while charging the battery
4. When the battery is full, the charge indicator light shows "GREEN". Please plug out the charger.

PRECAUTION:

1. Do not use the nebulizer while charging. The nebulizer cannot operate while charging.
2. To ensure safe using, the adapter and electronic equipment with USB port should meet the electrical safety standard IEC 60601-1.
3. The device is a BF type applied parts medical device.
4. When the battery is full, the charge indicator light shows "GREEN". Please plug out the charger.
5. The device did not use defibrillation-protected applied part.

Battery Recharge


6. The device is not intended for use in oxygen rich environment.
7. The device is a continuous operating device at steady flow rate, no flow rate adjustment available, requires manually turn on and off.
8. Always unplug the device from the power outlet after use and before cleaning.
9. Do not misuse the power cord. Do not wind the power cord around the compressor Always use the cable line with qualified power adapter.
10. Plug the device into the appropriate voltage outlet. Do not overload power outlets or use extension cords.
11. Do not operate the device with a damaged USB port or damaged cable line, may experience electrical leakage/shock.
12. Keep the charging cable line away from heated surfaces.
13. Do not plug or unplug the AC adapter into the electrical outlet with wet hands.
14. In order to prevent strangulation by the charging cable line, or injuries caused by mishandling, please pay close supervision when this device is being used by, on or near children or invalids.
15. Make sure the cables line use to charge is working fine, if not stop using the device. Contact local distributor for exchange or refund service.
16. Do not overload power outlets. Plug the AC adapter into the appropriate voltage outlet.
17. Do not use extension cords. Plug the AC adapter directly into the electrical outlet. Do not pull the power cord of AC adapter strongly.

Note:

AC adapter can be worked at 100 - 240V ~ 350mA, 50/60 Hz.

Change Air Filter

1. Pull the air filter cover off the device.
2. Change the air filter.
3. Put the air filter cover back in place.

CAUTION:

1. Replace the nebulizer after a long period of inactivity, in cases where it shows deformities, breakage, or when the atomizer head is obstructed by dry medicine, dust, etc.
2. If the filter shows a grey or brown colour or is wet, replace it. Extract the filter and replace it with a new one. Do not try to clean the filter for reusing it. The air filter shall not be serviced or maintained while in use with a patient. Only use original filters! Do not use the device without filter!
3. Replace the accessories/air filter routinely. Fail to do so, may derive bacterial infection or degradation of the performance.
4. Make sure that the air filter is clean. If the air filter has been used for more than 90 days, replace it with a new one.
5. Do not block the air filter cover.

Cleaning and Disinfecting

- Thoroughly clean all components to remove medication residuals and possible impurities after each treatment.
- Use a soft and dry cloth with non-abrasive cleaners to clean the unit.
- Make sure that the internal parts of the device are not in contact with liquids and that the charging cable is disconnected.
- Clean the medication tank, vaporized head and inhalation top after each use. Disinfect the accessories daily.
- Do not immerse the nebulizer main unit in water or other liquid.
- Do not distort the vaporizer head, or jab the nozzle of the medication tank with a pin, brush or any sharp object.
- Do not use a microwave oven, dish dryer or hair dryer to dry the device or the parts.
- Do not use an autoclave, EOG gas sterilizer or low temperature plasma sterilizer.
- When disinfecting parts by boiling, make sure that the container does not boil dry. Otherwise it may result in a fire and damage the parts.
- Do not boil vaporizer head together with objects other than applicable nebulizer accessories.

Cleaning

1.Cleaning Nebulizer Kit (Inhalation top, medication tank), Mouthpiece and Vaporizer head

- a. Disassemble all components except the main unit.
- b. Prepare mild detergent.
- c. Rinse all components with regular cold drinking water (around 20°C) for approximately two minutes.
- d. Use mild detergent to clean each part. The cleaning time is about five minutes.
- e. After cleaning with mild detergent, rinse with cold drinking water (about 20°C) for about two minutes.
- f. After the cleaning operation is completed, place all parts in a room temperature environment to completely air dry for about two hours.

Cleaning and Disinfecting

Mouthpiece is intended for daily cleaning if used and is valid for three months and a cleaning limitation of 1000 cycles.

2.Nebulizer Device

a. Wipe the device with a clean towel moistened with drinking quality water.

b. Wipe the device dry with a clean, dry towel.

Recommend to wipe the device on a daily bases.

Disinfecting

1.Disinfecting Nebulizer Kit (Inhalation top, medication tank), Mouthpiece and Vaporizer head

a. Prepare common household cooking pot.

b. Boil tap water and allow it to cool naturally to 70°C or below.

c. Put the Mouthpiece/Inhalation Top/Medication Tank/Vaporizer Head into the aforementioned warm water.

d. Soak for 5 minutes.

e. When the time is up, remove the accessories from the water and let them air dry completely at room temperature for about two hours.

Mouthpiece is intended for daily disinfecting if used and is valid for three months and a cleaning limitation of 1000 cycles.

2.Mask

a. Prepare alcohol swabs.

b. Use an alcohol swab to wipe the face shield, covering the entire surface area of the face shield.

Mask is intended for daily cleaning if used and is valid for three months and a disinfecting limitation of 100 cycles.

CAUTION:

Do not wash or clean the air filter. If the air filter becomes wet, replace it. Damp air filters can cause blockages.

Technical Specification

Thoroughly clean all components to remove medication residuals

Type: Compressor nebulizer

Particle Size: 2 - 5 μm

Nebulization Rate: approx. 0.25 + 0.15/-0.05 ml/ min (0.2~0.4 ml/min)

Noise: $\leq 65 \text{ dB(A)}$

Medication Cup Capacity°C max. 5ml

Power Supply: Rechargeable Li-ion Polymer Battery.

Enclosure Degree of Ingress Protection: IP22

Dimension: 195(L) x 72.5(W) x 54(H) mm

Weight: 163g $\pm 5\%$

Operating conditions:

+10 ° C to +40 ° C (50 - 104 °F)

10% to 95% relative humidity (non-condensing)

Atmospheric pressure: 700-1060 hPa

Storage conditions:

-20 ° C to +60 ° C (-4 - 140 °F)

10% to 95% relative humidity (non-condensing)

Atmospheric pressure: 700-1060 hPa

Battery Operation Time: estimated 1 hr

Lifetime: 400 hrs

Shelf Life: 3 years

Content: Device x 1, Manual x 1, USB Type-C cable,
Adult Mask x 1, Child Mask x 1, Mouthpiece x 1, FILTER x 5

Consult your physician or doctor if

You have the following medical symptoms or medical history

What may interact with this medication? Albuterol or Salbutamol

- Amphetamines (e.g., dextroamphetamine, lisdexamfetamine)
- Antiarrhythmics (e.g., amiodarone, dronedarone, procainamide)
- Antipsychotics (e.g., chlorpromazine, haloperidol, ziprasidone)
- Atomoxetine
- Beta-blockers (e.g., carvedilol, propranolol, labetalol, nadolol, sotalol)
- Caffeine
- Cannabis
- Chloroquine
- Cisapride
- Decongestant eye drops and nose sprays (e.g., naphazoline, oxymetazoline, xylometazoline)
- Disopyramide
- Diuretics
- Epinephrine
- Fast-acting bronchodilators (e.g., salbutamol, terbutaline)
- Linezolid
- Long-acting bronchodilators (e.g., formoterol, salmeterol)
- Medications for colds
- Medications for depression or for emotional or psychotic conditions
- Medications for weight loss including some herbal products
- Methadone
- Methylphenidate
- Midodrine
- Modafinil
- Monoamine oxidase inhibitors (e.g., phenelzine, tranylcypromine, moclobemide, rasagiline, selegiline)
- Nabilone
- Ozanimod
- Pentamidine

Consult your physician or doctor if

- Protein kinase inhibitors (e.g., lenvatinib, selpercatinib, vandetanib)
- Quinine
- Some antibiotics like clarithromycin, erythromycin, levofloxacin, and linezolid
- Some heart medications
- Steroid hormones like dexamethasone, cortisone, hydrocortisone
- Theophylline
- Thyroid hormones
- Tricyclic antidepressants (e.g., amitriptyline, desipramine, nortriptyline)

What may interact with this medication? Terbutaline

Do not take the medication with any of the following:

- Adrenalin (epinephrine)
- amitriptyline
- Atarax (hydroxyzine)
- doxepin
- hydrochlorothiazide
- labetalol
- Lasix (furosemide)
- metoprolol
- monoamine oxidase inhibitors
- phenelzine
- Phenergan (promethazine)
- Pitocin (oxytocin)
- selegiline
- Symbicort (budesonide / formoterol)
- tricyclic antidepressants,
- Zofran (ondansetron)

If you experience these side effects

Contact your physician or doctor immediately; your pharmacist may be able to advise you on managing side effects.

- Change in taste sensation or a bad taste in the mouth
- Changes in blood pressure
- Common cold (runny nose, nasal congestion, sore throat, Coughing)
- Constipation
- Cramps weakness,
- Diarrhea
- Difficult or painful urination
- Dizziness
- Dry mouth
- Feeling weak,
- Symptoms of glaucoma (e.g., blurred vision, seeing halos of bright colours around lights, red eyes, increased pressure in your eyes, eye pain or discomfort)
- Headache
- Hoarseness
- Loss of appetite
- Muscle pain, weakness, or spasms
- Mental disorder;
- Nausea
- Nervousness
- Respiratory infection
- Sweating
- Thrush—white patches in the mouth
- Tremor
- Trouble in swallowing or breathing
- Trouble sleeping
- Voice changes
- Vomiting
- Wheezing after inhalation;

Seek immediate medical attention

If any of the following occur, stop taking the medication:

- sudden blurred vision or eye pain
- fast or irregular heartbeat
- increased shortness of breath or wheezing, breathing difficulties or shortness of breath becomes worse
- signs chest pain or heart attack (e.g., sudden chest pain or pain radiating to back, down arm, jaw; sensation of fullness of the chest; nausea; vomiting; sweating; anxiety)
- signs of a serious allergic reaction (e.g., abdominal cramps, difficulty breathing, nausea and vomiting, or swelling of the face and throat)
- skin rash or hives
- increased Medication Dosage Required or you don't get the same benefit from your medicine as you did before;
- difficulty with urination

Troubleshooting

In case of an unexpected behavior of your device check with the troubleshooting table if you find the issue in the column “Problem” and follow the instructions in the column “Solution”.

Problem	Cause	Solution
After turning on the power switch, the device does not move.	The device is out of power.	Please charge the battery.
	The AC adapter is not plugged into an electrical outlet.	Turn the power switch off. Plug the AC adapter into an electrical outlet. Turn the device on.
No nebulization or low nebulization rate when the power is on.	There is no medicine in the medicine tank.	Add the correct amount of prescription medication to medicine tank.
	Improper amount of medicine in the medicine tank.	Make sure the nebulizer kit is correctly assembled and the inhalation accessory is correctly attached.
	The nebulizer kit is not correctly assembled.	Clean and disinfect the nebulizer kit to remove the blockage.
	The vaporizer head is blocked.	Hold the nebulizer kit correctly. Do not tilt the nebulizer kit.
	The nebulizer kit is tilted at an incorrect angle.	Replace the air filter with a new clean air filter.

EMC Tables

NE10 is intended for use in the electromagnetic environment specified below. The customer or the user of **NE10** must make sure that it is used in such an environment.

Guidance and manufacturer's declaration - Electromagnetic emissions

Phenomenon	Professional healthcare facility environment a)	HOME HEALTHCARE ENVIRONMENT a)
Conducted and radiated RF MISSIONS	CISPR 11 Group 1 Class A (Not BLE Function) Group 2 Class A (With BLE Function)	CISPR 11 Group 1 Class B (Not BLE Function) Group 2 Class B (With BLE Function)
Harmonic distortion	Not applicable b) (Note: Power by Battery or DC Input) Only the AC input needs to be tested	
Voltage fluctuations and flickering	Not applicable b) (Note: Power by Battery or DC Input) Only the AC input needs to be tested	

a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive acceptance limits of Group 1 Class B (CISPR 11) have been considered and applied. The equipment is suitable for use in the mentioned environments when directly connected to the Public Mains Network.

b) The test is not applicable in this environment unless the ME EQUIPMENT and ME SYSTEM used will be connected to the PUBLIC MAINS NETWORK and the power input is otherwise within the scope of the Basic EMC standard.

Guidance and manufacturer's declaration - Electromagnetic immunity - Enclosure port

Phenomenon	Basic EMC standard or test method	Immunity test levels	
		Professional healthcare facility environment	HOME HEALTHCARE

EMC Tables

ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8kV contact ± 2 kV, ±4kV ±, ±8 kV, ±15 kV air	
Radiated RF EM fields	IEC 61000-4-3	a)	10 V/m b) 80MHz - 2.7 GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	COMPLIANT NOTE: Further information about distances to be maintained between portable and mobile RF communications equipment (transmitters) and the NE10 can be requested from supplier using the contact information provided in this manual. However, it is advisable to keep the equipment at an adequate distance of, at least, 0.5 m from mobile phones or other RF communications transmitters to minimize possible interference.	
RATED power frequency magnetic fields.	IEC 61000-4-8	30 A/m c) 50 Hz or 60 Hz	
a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive IMMUNITY acceptance limits have been considered and applied. b) Before modulation is applied. c) This test level assumes a minimum distance of at least 15 cm between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic fields.			

Explanation of Symbols

	The CE marking with the Registration Number of the Notified Body. This denotes the compliance of Regulation (EU) 2017/745
	Medical Device
	Manufacturer
	Authorized representative in the European Community
	Date of manufacture (YYYY-MM-DD or YYYY-MM)
	Batch code (YYMMWWWW)
	Serial number (YYMWWWWXXXXX)
	Keep dry
	Indicates a medical device that needs protection from light sources.
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation
	Caution
	Consult the instruction for use
	Disposal information: Should you wish to dispose of the article, do so in accordance with current regulations. Details are available from your local authority. WEEE 2012/19/EU Directives

Explanation of Symbols

RoHS	This product fulfilling the requirements of the RoHS Directive 2011/65/EU.
REACH	This product fulfilling the requirements of the REACH Directive EC 1907/2006 and its amendments, do not contain Substances of Very High Concern in concentration above the limit of 0.1 %. No substance(s) is/are present in the parts of the product above the concentration of 0.1 % weight by weight.
	"ON/OFF" (push-push)
	Device classification type BF
IP22	This product meets the basic safety and essential performance requirements indicated in the IP22 conditioning test (protection against solid foreign objects of 12.5mm Ø and greater and against vertically falling water drops when enclosure tilted up to 15°)
	The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.
	Importer
	Distributor
	Model Number
	Catalogue Number
	Country of Manufacturer
	Unique Device Identifier
	WARNING Projectile Hazard The Nebulizer is MR Unsafe.

CONDITIONS FOR SERVICE UNDER GUARANTEE

1. SANICO Srl guarantees the absence of faults and mechanical defects and the completeness of the product when it is sold.
2. Repairs can be done only at Service Centres authorized by SANICO Srl or by the Authorized Distributor in the country where you have purchased the appliance.
3. The duration of the period of guarantee on the materials and on manufacturing and operating defects is 2 years with effect from the date of sale, subject to all the conditions for the Purchaser shown in the Instructions for use. The obligations of guarantee do not cover expendable materials: electrode cups, elastic bands, power supply cables etc. nor regular checks and maintenance (calibration).
4. Complaints relative to the product must be presented directly to the authorized service centre. The service centre, on the basis of an agreement with the Vendor, is competent to establish the cause and nature of the fault. The product will be replaced only if it cannot be repaired.
5. The owner will provide to deliver the product to the after-sales service center, as per the delivery instructions that can be requested by e-mail: info@sanicare.it.
6. The time necessary for the repair under guarantee will be subject to the degree of seriousness of the fault.
7. The responsibility of SANICO Srl under these conditions of guarantee will be limited to the amount paid by the Purchaser for the product in question.

The guarantee on the product will not be applied in the following cases

- certificate of guarantee filled in incorrectly or falsified; serial number (if applicable for the specific product) damaged or removed;
- absence of commercial and financial documents (cash receipt, delivery note) proving the purchase;
- opening or repair of the product by the Purchaser or by an unauthorized Service Centre;
- defects following failure by the Purchaser to observe the rules for use,

- storage and transport of the product shown in the relative technical and regulations documentation, in which case the regulations of the aforementioned documentation which have been infringed must be specified;
- penetration of foreign bodies or liquids inside the product;
- presence on the product of internal or external mechanical defects (cracks, signs of a fall, chipped parts, etc.) consequent to the incorrect use, installation or transport of the product;
- components or accessories subject to rapid wear and tear (electrode cups, elastic bands);
- defects caused by external agents (dirt visible externally and internally, traces of corrosion);
- defects due to causes independent of the manufacturer, such as power supply voltage step changes, natural phenomena or disasters, fires, actions by animals, insects etc.
- complaints for defects found by the Purchaser that are made after the expiry of the guarantee period or any other period established by the law for reporting any defects;
- use of the product by the Purchaser not for personal (domestic) use but for commercial and profit-making purposes.

Conditions for the return of a defective product

Defective products will be accepted only if clean and complete, in the original packaging, together with the certificate of guarantee correctly filled in and the commercial and financial documents (cash receipt of the store, delivery note) that prove the purchase. Please attach a description of the problem found to the appliance.

After Sales Service Center, please contact:

e-mail: info@sanicare.it

www.sanicare.it

Before sending the appliance, please contact us by phone or e-mail to agree on how the product is to be returned

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(se disponibile/if available)

Data d'acquisto/Purchase date: _____

Apparecchio in Garanzia/Under warranty:

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WARRANTY IS VALID ONLY IF ACCOMPANIED BY INVOICE/TICKET.

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DISPOSAL OF YOUR OLD APPLIANCE

When this crossed-out wheeled bin symbol is attached to a product, it means the product is covered by the European Directive 2002/96/EC. All electrical and electronic products should be disposed off separately from the municipal waste stream via designated collection facilities appointed by the government or the local authorities. The correct disposal of your old appliance will help prevent potential negative consequences for the environment and human health. For more detailed information about disposal of your old appliance, please contact your city office, waste disposal service or the shop where you purchased the product.

AVITA Corporation
9F., No. 78, Sec. 1, Kwang Fu Rd., San Chung Dist.,
New Taipei City 24158, Taiwan, China

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MDSS GmbH